


MITTATAULUKKO

		Koko					
Mittauspisteet		I	II	III	IV	V	VI
Ympärysmitta (cm)	cG silikonilla, normaali reiden ympärysm	43-48	46-52	49-56	53-60	56-64	60-68
	cD silikonilla, suuri ympärysm	49-57	53-62	57-67	61-72	65-77	69-81
	cG ilman silikonia	43-57	46-62	49-67	53-72	56-77	60-81
	cF	39-52	41-56	44-60	47-65	50-69	53-73
	cE	30-37	33-40	36-43	37-45	39-48	41-51
	cD	27-33	29-36	32-39	34-42	36-45	38-46
	cC	28-34	30-37	33-40	35-43	37-46	39-49
	cB1	23-27	24-29	26-32	29-35	31-37	33-39
	cB	18-20	20-22	22-24	24-26	26-28	28-30
	cY	26-31	26-33	29-35	31-37	32-38	33-40
	cA	17-22	19-24	21-26	23-29	25-32	27-34

Tehdaskissa on 5cm leveä silikonit.



Pitusmitta:

AD
 ID lyhyt: 34 - 50 cm
 ID normaali: 39 - 45 cm

AG
 IG lyhyt: 60 - 69 cm
 IG normaali: 70 - 80 cm
 IG pitkä: 81 - 90 cm (kokoaa
 JORST® Classic)

AT/AG-T
 IK lyhyt: 65 - 75 cm
 IK normaali: 75 - 85 cm